



#### ABDOMINAL WOUND

1. Abdominal bandage with transparent plastic foil.
2. Place the abdominal bandage firmly on the wound.
3. Wrap the elastic bandage around the body...
4. To cover up all edges of the wound pad.
5. Wrap the elastic bandage firmly...
6. Secure the bandage with the closure bar.

#### AMPUTATIONS

1. Wrap the wound pad around the limb.
2. Wrap the elastic bandage around the wound pad.
3. Twist the elastic bandage a few turns to create a cord...
4. And wrap the cord around the limb...
5. Insert the closure bar between the cord and rotate, to create a tourniquet effect...
6. Secure the bandage with the closure bar.